

TABELA DE COMPATIBILIDADE
MICROPLACAS PARA PCR X TERMOCICLADORES



		PCR-96-FS-C	PCR-96M2-HS-C	PCR-96-MB-C	PCR-96-AB-C	PCR-96-C	PCR-48-C	PCR-24-C
Axygen	Maxygene - THERM-1000	●	●	●	●	●	●	●
AmpliTherm	96 poços	●	●	●	●	●	●	●
	96 poços com gradiente	●	●	●	●	●	●	●
Applied Biosystems	Veriti	●	●	●	●	●	●	●
	9600	●	●	●	●	●	●	●
	9700	●	●	●	●	●	●	●
	9800	●	●	●	●	●	●	●
	2720	●	●	●	●	●	●	●
Biocycler	MJ96+	●	●	●	●	●	●	●
	MJ96G	●	●	●	●	●	●	●
Biometra	Uno	●	●	●	●	●	●	●
	Uno II	●	●	●	●	●	●	●
	T1 Thermalcycler	●	●	●	●	●	●	●
	Tgradient	●	●	●	●	●	●	●
	Trobot	●	●	●	●	●	●	●
BioRad	ICycler	●	●	●	●	●	●	●
	MyCycler	●	●	●	●	●	●	●
Corbett Research	PalmCycler 96	●	●	●	●	●	●	●
	PalmCycler 384	●	●	●	●	●	●	●
Eppendorf	Mastercycler	●	●	●	●	●	●	●
	Mastercycler ep	●	●	●	●	●	●	●
	Mastercycler M384	●	●	●	●	●	●	●
Ericom	Single Block	●	●	●	●	●	●	●
	TwinBlock	●	●	●	●	●	●	●
	Deltacycler I	●	●	●	●	●	●	●
Flexi	Gene	●	●	●	●	●	●	●
	Genius	●	●	●	●	●	●	●
MJ Research	PTC-100	●	●	●	●	●	●	●
	PTC-200	●	●	●	●	●	●	●
	PTC-225 Tetrad	●	●	●	●	●	●	●
	Dyad/Disciple	●	●	●	●	●	●	●
MWG	Primus 96	●	●	●	●	●	●	●
	Primus 384	●	●	●	●	●	●	●
Stratagene	Robocycler	●	●	●	●	●	●	
TaKaRa	TP240	●	●	●	●	●	●	●
	TP3000	●	●	●	●	●	●	●
Techne	Touchgene Gradient (TC512)	●	●	●	●	●	●	●
	Genius (TC412)	●	●	●	●	●	●	●
	Touchgene X	●	●	●	●	●	●	●
ThermoHybaid	PCR Sprint	●	●	●	●	●	●	●
	PCR Express	●	●	●	●	●	●	●
	Multiblock System	●	●	●	●	●	●	●
	Touchdown	●	●	●	●	●	●	●
	Omnigene	●	●	●	●	●	●	●
	Omn-E	●	●	●	●	●	●	●
	PxE	●	●	●	●	●	●	●
	Px2	●	●	●	●	●	●	●

TABELA DE COMPATIBILIDADE MICROPLACAS PARA PCR X TERMOCICLADORES DE TEMPO REAL



		PCR-96-FS-C	PCR-96M2-HS-C	PCR-96-MB-C	PCR-96-AB-C	PCR-96-C	PCR-48-C	PCR-24-C
Applied Biosystems	5700	●	●	●	●	●	●	●
	7000	●	●	●	●	●	●	●
	7300	●	●	●	●	●	●	●
	7500	●	●	●	●	●	●	●
	7700	●	●	●	●	●	●	●
	7900HT	●	●	●	●	●	●	●
BioRad	ICycler	●	●	●	●	●	●	●
	MyiQ	●	●	●	●	●	●	●
	IQ5	●	●	●	●	●	●	●
Eppendorf	MasterCycler ep Realplex	●	●	●	●	●	●	●
MJ Research	Opticon	●	●	●	●	●	●	●
	Opticon 2	●	●	●	●	●	●	●
	Chromo 4	●	●	●	●	●	●	●
Stratagene	Mx4000	●	●	●	●	●	●	●
	Mx3000	●	●	●	●	●	●	●
Techne	Quantica	●	●	●	●	●	●	●

TABELA DE COMPATIBILIDADE MICROPLACAS PARA PCR X SEQUENCIADORES

		PCR-96-FS-C	PCR-96M2-HS-C	PCR-96-MB-C	PCR-96-AB-C	PCR-96-C	PCR-48-C	PCR-24-C
Amersham	MegaBACE 500	●	●	●	●	●	●	●
	MegaBACE 1000	●	●	●	●	●	●	●
	MegaBACE 4000	●	●	●	●	●	●	●
Applied Biosystems	3100	●	●	●	●	●	●	●
	3130	●	●	●	●	●	●	●
	3700	●	●	●	●	●	●	●
	3730/3730x	●	●	●	●	●	●	●
Beckman	CEQ	●	●	●	●	●	●	●
MJ Research	BaseStation	●	●	●	●	●	●	●
Transgenomic	WAVE System	●	●	●	●	●	●	●

Distribuidor